

Laboratory and Operating Rooms Intercom Systems

Model IC-FA

- It was specifically designed for laboratories and companies that have a need to communicate with cleanroom, white areas and / or sectors where operators pack products and they cannot press any button for health or safety requirement.
- Front built in stainless steel with rounded edges buttons to prevent accumulation of dust.
- The front, electronic board, speaker and high sensitivity omnidirectional microphone are mounted on the wall or in a plastic wall box . Roker line Model PR 1000/81.
Dimensions: 115 x 158 x 60
Stainless steel front: 125x170 mm.

The model IC F. A. semi-handsfree allows:

- Communication from the hall IC model CA with the sterile area, just by pressing one button. From the sterile area they respond "hands-free".
- If they require to call from the sterile area IC model RTA /S to the hallway, they press the button only once and the light will activate intermittently.



Laboratory and Operating Rooms Intercom Systems

Model IC-RTA/S



Model IC-CA



Luminous device

Model IC-134/RB



Red intermittently light + buzzer

www.intercron.com

Telephone number / Fax: 54 11 4545-5400 (telephone line rollover) - E-mail: info@intercron.com
Address: Valdenegro 3201 - (C1431GAM) - Capital Federal - Argentina